

Active II/360 Water Birth Pool


Originally designed with the help of a health service ergonomist and midwives in 1995 this large, luxurious birthing pool will help you create a world-class facility for labour and birth in water.

We've transformed the Active II into a freestanding pool and added a specially designed Water Column to maximise hygiene, water safety and provide support for mothers and midwives.



Specifications

Length:	1965mm
Width:	1650mm
Height:	770mm
Depth:	640mm
Weight (empty):	100kg
Capacity:	560 litres
Material:	Ficore® composite
Guarantee:	Lifetime



Active II/360 Water Birth Pool

ACTIVE
BIRTHPOOLS

The bull-nose shaped rim provides firm hand-hold for support in upright positions

Multi-colour LED lighting and Bluetooth sound system allow mothers to personalise the ambiance of the room

The sturdy column provides support for upright birth positions and a plinth for midwives to lean on

The keyhole shape channel in the back rim protects mothers from hitting their head on the spout

The recessed hand grips provide the ideal support for upright positions such as squatting

The wide rim, grab rail and step unit make it simple and safe for mothers to get in and out of the pool

Designed to safeguard mothers and midwives from risk of injury and cross infection

The rounded, bull nose-shaped rim provides dynamic support for upright positions

Specialised material makes our pools less slippery, tactile and warm to the touch

This large, luxurious pool gives mothers and midwives complete freedom of movement both in and around the pool

The extra hard, dense material is highly resistant to damage, bacteria and disinfectants

Midwives sit naturally, in comfort with their legs under the pool and arms resting on the extra-wide rim

The concave, cutaway surround enables portable hoists to be utilised for emergency evacuations

The innovative Safety Seat provides support for mothers and gives midwives safe, practical options for monitoring, and emergency evacuations

Mothers have optimum space, depth and support to gain maximum benefit from water immersion

The unique Labour Support Seat is tilted forward to encourage the pelvis to spread and the birth canal to open.

