

Venus 360 Water Birth Pool

The mid-size Venus 360 is the perfect choice for rooms where end users want to have the pool away from the wall.


It combines intelligent design with superior materials and build quality to create a birthing pool that delivers superior value, outstanding performance, and unparalleled safety.

We've increased access around the back of the pool and developed the Water Column to create a birth pool that is safer and more user-friendly.



Specifications

Length:	2065mm
Width:	1350mm
Height:	770mm
Depth:	640mm
Weight (empty):	90kg
Capacity:	510 litres
Material:	Ficore® composite
Guarantee:	Lifetime



Venus 360 Water Birth Pool

ACTIVE
BIRTHPOOLS

Multi-colour lighting and Blue-tooth sound enable mothers to personalise the ambience of the room

The Water Column separates the plumbing from the pool to maximise hygiene and water safety

The 'keyhole' shaped channel protects the mother from hitting her head on the spout

The recessed hand grips are bonded into fabric of the pool to deny bacteria to a place to grow

The extra-wide rim, grab rail and step unit make it easy for mothers get in and out of the pool

The combination of double wall construction and ultra-dense material reduces heat loss to .7 degrees per hour

Provides freedom of movement for mothers and midwives to move naturally both in and around the pool

Seamless one-piece construction in Ficare composite makes the pool more durable and easier to clean

The Labour support Seat encourages optimal positioning and gives mothers a place to rest and welcome their baby

Mothers use the extra-wide rim to rest on and to safely enter and leave the pool

Easier to clean and maintain and can be easily repaired on site like new if damaged

The concave, cut-away surround makes it easier and more comfortable for midwives to be close to the mother

Optimal space and depth for mothers to move and find most beneficial positions

The bull-nose shaped rim provides mothers with support in the wide range of postures natural to labour and birth

