


# TECHNICAL SPECIFICATIONS

**ACTIVE**  
BIRTHPOOLS  
www.activebirthpools.com

## Princess II/360 Birth Pool

<b>Material</b>	Ficore composite
<b>Capacity</b>	425 liters / 112 Gallons (4" below rim)
<b>Height</b>	770mm / 30"
<b>Width</b>	1100mm / 43"
<b>Length</b>	1900mm / 75"
<b>Water Depth</b>	640mm / 25" (4" below rim)
<b>Weight (empty)</b>	70kg / 154 lbs
<b>Shipping Crate</b>	167 x 117 x 89cm / 66" x 46" x 35"
<b>Shipping Weight</b>	145 kg / 320 lbs
<b>Birth Stool</b>	1 - Submerged
<b>Handles</b>	None
<b>Drain Mechanism</b>	40mm / 1 1/2" Mechanical Pop-up—Brass with cable mechanism, rim mounted control, removable and height adjustable chrome stopper
<b>Waste Plumbing</b>	40mm / 1 1/2" PVC - Shallow P-trap included - thru-wall (or thru-floor connection - omit supplied p-trap)
<b>Access Panel</b>	2 - right and left rear of pool - 10mm hex key lock
<b>Chromotherapy Light</b>	1 - 9 light, 11 color LED: White, light blue, dark blue, purple, magenta, red, pink, orange, yellow, light green, dark green
<b>Transformer</b>	110VAC to 12VDC 0.5A - CSA and UL Certified - Nema 5-15 Plug
<b>Light Controls</b>	Rim mounted - chrome, sealed, four button - Mounted on right side unless specified
<b>Step Stool</b>	18" x 12" x 6", 8.5 lbs, Ficore composite

www.activebirthpools.com



NOTE: Dimensions are in millimeters if not indicated otherwise.

# Princess II/360 Water Birth Pool

ACTIVE  
BIRTHPOOLS

The Labour Support Seat encourages the pelvis to open and provides support in upright postures

The Water Column separates the plumbing from the pool to guarantee highest level of hygiene and water safety

The keyhole shaped channel in the rim protects mothers from hitting their head on the spout

Integral LED Lighting and Bluetooth Sound allow mothers to personalise the ambiance of the room

The step unit and grab rail help make it easy for mothers to safely get in and out of the pool

The wrap-around rim provides support for mothers in the full range of positions natural to labour and birth

The highly polished finish is less slippery, super smooth, tactile and warm to the touch

360 degrees of movement around and in the pool for mothers and midwives

The Lumbar Support provides comfort and support for mothers in upright postures like squatting

High density and superior strength minimise heat loss and make our pools more durable and economical

Though compact, the Princess Pool is surprisingly spacious for mothers to move freely in comfort

Designed to safely facilitate emergency evacuations with or without a portable hoist

Easy to clean, resistant to bacteria and repairable on site, like new if damaged

There's plenty of space under the rim for midwives to get close to the mother in comfort

